

訂 正

以下の論文に一部誤りがありましたので訂正いたします。

感染症学雑誌 第 89 巻 第 5 号 574~578 頁

タイトル：新規に開発された respiratory syncytial virus 迅速診断キットの有用性評価

著者：黒木 春郎 他

574 ページ

要旨 13 行目

(誤)：A 型 (Wild 株) 3.0×10^2 TCID₅₀/mL, B 型 (Long 株)

↓

(正)：A 型 (Long 株) 3.0×10^2 TCID₅₀/mL, B 型 (Wild 株)

576 ページ

左側 2 行目

(誤)：100.0% (71/71), 94.2% (131/139)

↓

(正)：89.9% (71/79), 100.0% (131/131)

左側 5 行目

(誤)：90.1% (64/71), 89.2% (124/139)

↓

(正)：81.0% (64/79), 94.7% (124/131)

左側 9 行目

(誤)：100.0% (50/50), 86.1% (93/108)

↓

(正)：76.9% (50/65), 100.0% (93/93)

左側 19 行目

(誤)：95.8% (46/48), 96.5% (83/86)

↓

(正)：93.9% (46/49), 97.6% (83/85)

右側 7 行目

(誤)：97.9% (47/48), 97.6% (83/85)

↓

(正)：95.9% (47/49), 98.8% (83/84)

Table 1

(誤)：

Nasal swab	RT-PCR			Concordance Rate	
	Positive	Negative	Total		
IC-A	Positive	71	0	71	Positive: 100.0%
	Negative	8	131	139	Negative: 94.2%
	Total	79	131	210	Total: 96.2%
IC-B	Positive	64	7	71	Positive: 90.1%
	Negative	15	124	139	Negative: 89.2%
	Total	79	131	210	Total: 89.5%
IC-C	Positive	50	0	50	Positive: 100.0%
	Negative	15	93	108	Negative: 86.1%
	Total	65	93	158	Total: 90.5%

↓
(正)：

Nasal swab	RT-PCR			Concordance Rate	
	Positive	Negative	Total		
IC-A	Positive	71	0	71	Positive: 89.9%
	Negative	8	131	139	Negative: 100.0%
	Total	79	131	210	Total: 96.2%
IC-B	Positive	64	7	71	Positive: 81.0%
	Negative	15	124	139	Negative: 94.7%
	Total	79	131	210	Total: 89.5%
IC-C	Positive	50	0	50	Positive: 76.9%
	Negative	15	93	108	Negative: 100.0%
	Total	65	93	158	Total: 90.5%

Table 2

(誤)：

Nasal swab	RT-PCR			Concordance Rate	
	Positive	Negative	Total		
IC-A	Positive	46	2	48	Positive: 95.8%
	Negative	3	83	86	Negative: 96.5%
	Total	49	85	134	Total: 96.3%
IC-B	Positive	45	4	49	Positive: 91.8%
	Negative	4	80	84	Negative: 95.2%
	Total	49	84	133	Total: 94.0%
IC-C	Positive	47	1	48	Positive: 97.9%
	Negative	2	83	85	Negative: 97.6%
	Total	49	84	133	Total: 97.7%

↓
(正)：

Nasopharyngeal aspirate	RT-PCR			Concordance Rate	
	Positive	Negative	Total		
IC-A	Positive	46	2	48	Positive: 93.9%
	Negative	3	83	86	Negative: 97.6%
	Total	49	85	134	Total: 96.3%
IC-B	Positive	45	4	49	Positive: 91.8%
	Negative	4	80	84	Negative: 95.2%
	Total	49	84	133	Total: 94.0%
IC-C	Positive	47	1	48	Positive: 95.9%
	Negative	2	83	85	Negative: 98.8%
	Total	49	84	133	Total: 97.7%